



# AURELIA® VINTAGE

Powdered Latex Disposable Gloves

## AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND  
sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building 1000,  
Gateway Business Park,  
New Mallow Road, Cork,  
Ireland

[aureliagloves.ie](http://aureliagloves.ie)

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

 10 x 100



















XS	S	M	L	XL
28825	28826	28827	28828	28829

Sizes & Dimensions	Single wall thickness	Physical Properties
--------------------	-----------------------	---------------------

 XS 70+10 mm S 80+10 mm M 95+10 mm L 110+10 mm XL > 110 mm	 min. 240 mm min. 240 mm min. 240 mm min. 240 mm min. 240 mm	 ① min. 0.08 mm ② min. 0.08 mm ③ min. 0.06 mm
--	--	---

	Unaged	Aged
Tensile Strength	Min. 18.0 MPa	Min. 14.0 MPa
Elongation	Min. 650%	Min. 500%
Strength at Break	Min. 6N	Min. 6N

 Latex	 Powdered	 3 Years
 Smooth	 5.5g +/- 0.2g	 Natural

