



AURELIA® ROBUST 9.0

5mil, Nitrile Powder Free Examination Gloves

| | |
|----------------------|--|
| Regulations | PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC |
| Quality Standards | ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems |
| Harmonised Standards | BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs. |

AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building 1000,
Gateway Business Park,
New Mallow Road, Cork,
Ireland

aureliagloves.ie



10 x 100



UK
CA
0120




XS S M L XL

96895 96896 96897 96898 96899

Sizes & Dimensions

Single wall thickness

Physical Properties



XS 70+10 mm
S 80+10 mm
M 95+10 mm
L 110+10 mm
XL > 110 mm



min. 240 mm
min. 240 mm
min. 240 mm
min. 240 mm
min. 240 mm



① min. 0.13 mm
② min. 0.10 mm
③ min. 0.08 mm

| | Unaged | Aged |
|-------------------|---------------|---------------|
| Tensile Strength | Min. 14.0 MPa | Min. 14.0 MPa |
| Elongation | Min. 500% | Min. 400% |
| Strength at Break | Min. 9N | Min. 9N |



Nitrile



Powder Free



5 Years



Finger Textured 5.5g +/- 0.2g



Ice Blue

