



AURELIA[®] ROBUST

4.5mil, Nitrile Powder Free Examination Gloves

| | |
|----------------------|--|
| Regulations | PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC |
| Quality Standards | ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems |
| Harmonised Standards | BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1184: 2002. Materials and articles in contact with foodstuffs. |

AURELIA[®] GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building 1000,
Gateway Business Park,
New Mallow Road, Cork,
Ireland

aureliagloves.ie



10 x 100

MD

UK
CA
0120

CE
0598

AQL
1.5

EN ISO 374-5: 2016
VIRUS

EN ISO 374-1: 2016/Type B
KPT



LATEX


XS S M L XL XXL

93895 93896 93897 93898 93899 93890

Sizes & Dimensions

Single wall thickness

Physical Properties



XS 70+10 mm
S 80+10 mm
M 95+10 mm
L 110+10 mm
XL > 110 mm



min. 240 mm
min. 240 mm
min. 240 mm
min. 240 mm
min. 240 mm



① min. 0.12 mm
② min. 0.09 mm
③ min. 0.05 mm

| | Unaged | Aged |
|-------------------|---------------|---------------|
| Tensile Strength | Min. 14.0 MPa | Min. 14.0 MPa |
| Elongation | Min. 500% | Min. 400% |
| Strength at Break | Min. 6N | Min. 6N |



Nitrile



Powder Free



5 Years



Micro



5.0g +/- 0.2g



Regular Blue

