



# AURELIA® SONIC 200

## 2.2mil, Nitrile Powder Free Examination Gloves

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1184: 2002. Materials and articles in contact with foodstuffs.

**AURELIA® GLOVES**  
A SUPERMAX HEALTHCARE BRAND  
sales@supermax.ie  
+353 15256724  
Units 1001 – 1012, Building 1000,  
Gateway Business Park,  
New Mallow Road, Cork,  
Ireland  
[aureliagloves.ie](http://aureliagloves.ie)

10 x 200

XS	S	M	L	XL
9377A5	9377A6	9377A7	9377A8	9377A9

Sizes & Dimensions	Single wall thickness	Physical Properties												
<p>XS 70+10 mm S 80+10 mm M 95+10 mm L 110+10 mm XL &gt; 110 mm</p> <p>min. 240 mm min. 240 mm min. 240 mm min. 240 mm min. 240 mm</p>	<p>① min. 0.09 mm ② min. 0.05 mm ③ min. 0.05 mm</p>	<table> <tr> <th></th><th>Unaged</th><th>Aged</th></tr> <tr> <td>Tensile Strength</td><td>Min. 14.0 MPa</td><td>Min. 14.0 MPa</td></tr> <tr> <td>Elongation</td><td>Min. 500%</td><td>Min. 400%</td></tr> <tr> <td>Strength at Break</td><td>Min. 6N</td><td>Min. 6N</td></tr> </table>		Unaged	Aged	Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa	Elongation	Min. 500%	Min. 400%	Strength at Break	Min. 6N	Min. 6N
	Unaged	Aged												
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa												
Elongation	Min. 500%	Min. 400%												
Strength at Break	Min. 6N	Min. 6N												

Nitrile
 Powder Free
 5 Years

Finger Textured
 3.0g +/- 0.2g
 Cobalt Blue

