




AURELIA® DELIGHT BLUE PD

Powdered Vinyl Examination Gloves



AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND
sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building 1000,
Gateway Business Park, New
Mallow Road, Cork,
Ireland

aureliagloves.ie

Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186:2002 Materials and articles in contact with foodstuffs.



10 x 100

S M L XL

38896 38897 38898 38899

Sizes & Dimensions

Single wall thickness

Physical Properties

	
S 80+10 mm	min. 240 mm
M 95+10 mm	min. 240 mm
L 110+10 mm	min. 240 mm
XL > 110 mm	min. 240 mm



- ① min. 0.05 mm
- ② min. 0.05 mm
- ③ min. 0.04 mm

	Unaged	Aged
Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa
Elongation	Min. 350%	Min. 350%
Strength at Break	Min. 3.6N	Min. 3.6N



Vinyl



Powdered



3 Years



Smooth



4.5g +/- 0.2g



Blue

