



# AURELIA® DELIGHT CLEAR PF

Premium Powder Free Vinyl Examination Gloves

**AURELIA® GLOVES**

A SUPERMAX HEALTHCARE BRAND

sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building 1000,  
Gateway Business Park,  
New Mallow Road, Cork,  
Ireland

[aureliagloves.ie](http://aureliagloves.ie)



|                      |                                                                                                                                                                                                         |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Regulations          | PPE Regulation (EU) 2016/425<br>Medical Device Regulation (EU) 2017/745, 93/42/EEC                                                                                                                      |
| Quality Standards    | ISO 9001:2015 Quality Management Systems<br>ISO 13485: 2016 Medical Devices Quality Management Systems                                                                                                  |
| Harmonised Standards | EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods.<br>BS EN 455 1-4 Medical gloves for single use.<br>BS EN 1186:2002 Materials and articles in contact with foodstuffs. |



10 x 100

**S M L XL**

38226 38227 38228 38229

## Sizes & Dimensions

## Single wall thickness

## Physical Properties

|    |           |             |
|----|-----------|-------------|
|    |           |             |
| S  | 80+10 mm  | min. 240 mm |
| M  | 95+10 mm  | min. 240 mm |
| L  | 110+10 mm | min. 240 mm |
| XL | > 110 mm  | min. 240 mm |



- ① min. 0.05 mm
- ② min. 0.05 mm
- ③ min. 0.04 mm

|                   | Unaged        | Aged          |
|-------------------|---------------|---------------|
| Tensile Strength  | Min. 11.0 MPa | Min. 11.0 MPa |
| Elongation        | Min. 350%     | Min. 350%     |
| Strength at Break | Min. 3.6N     | Min. 3.6N     |



Vinyl



Powder Free



3 Years



Smooth



4.5g +/- 0.2g



Clear

