



# AURELIA® VIBRANT

## Premium Natural Latex Disposable Gloves

### AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building  
1000, Gateway Business Park,  
New Mallow Road, Cork,  
Ireland

[aureliagloves.ie](http://aureliagloves.ie)

|                      |                                                                                                                                                                                                                                                                                                        |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Regulations          | PPE Regulation (EU) 2016/425<br>Medical Device Regulation (EU) 2017/745, 93/42/EEC                                                                                                                                                                                                                     |
| Quality Standards    | ISO 9001:2015 Quality Management Systems<br>ISO 13485: 2016 Medical Devices Quality Management Systems                                                                                                                                                                                                 |
| Harmonised Standards | BS EN ISO 374 - 1,2,4,5 Protective gloves against dangerous chemicals and microorganisms.<br>BS EN 455 1-4 Medical gloves for single use.<br>EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods.<br>BS EN 1186: 2002. Materials and articles in contact with foodstuffs. |



10 x 100



UK  
CA  
0120

CE  
0598



XS

S

M

L

XL

98265

98226

98227


98228

98229

#### Sizes & Dimensions

#### Single wall thickness

#### Physical Properties



XS 70+10 mm  
S 80+10 mm  
M 95+10 mm  
L 110+10 mm  
XL > 110 mm



min. 240 mm  
min. 240 mm  
min. 240 mm  
min. 240 mm  
min. 240 mm



① min. 0.13 mm  
② min. 0.10 mm  
③ min. 0.07 mm

|                   | Unaged        | Aged          |
|-------------------|---------------|---------------|
| Tensile Strength  | Min. 18.0 MPa | Min. 14.0 MPa |
| Elongation        | Min. 650%     | Min. 500%     |
| Strength at Break | Min. 6N       | Min. 6N       |



Latex



Powder Free



3 Years



Micro



5.7g +/- 0.2g



Natural

