



AURELIA® DELIGHT BLUE PF

Powder Free Blue Vinyl Disposable Gloves

AURELIA® DELIGHT BLUE PF

Powder Free Blue Vinyl Disposable Gloves




Vinyl



Powder Free



4.5g



Smooth Textured



Screen Compatible

 10 x 100

S	M	L	XL
38996	38997	38998	38999



BS EN ISO 374-5: 2016



VIRUS



BS EN ISO 374-1: 2016/Type B



KPT





Applications





Physical Properties

	Unaged	Aged
Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa
Elongation	Min. 350%	Min. 350%
Strength at Break	Min. 3.5N	Min. 3.5N

Regulations

Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

Innecase Dimensions:
235 x 125 x 60 mm

Outercase Dimensions:
325 x 255 x 250 m


EAN 13 Barcodes:
Inner
S: 955-5002-106126
M: 955-5002-106133
L: 955-5002-106140
XL: 955-5002-106154

Outer
S: 1955-5002-106123
M: 1955-5002-106130
L: 1955-5002-106147
XL: 1955-5002-106157

Features and Benefits


- Bright Blue** Ideal for food safety and preparation, the distinct colour helps easily identify contaminants and enhances visibility during handling.
- Food Contact Approved** Meets safety standards for direct food contact, ensuring hygienic handling and preventing contamination in food preparation environments.
- Smooth-Textured** Ensures a comfortable fit and effortless handling, offering a sleek surface for easy grip and reduced friction.

Single Wall Thickness



- 1 min. 0.07 mm
- 2 min. 0.05 mm
- 3 min. 0.03 mm

Sizes & Dimenensions



- | | |
|----|-----------|
| S | 80+10 mm |
| M | 95+10 mm |
| L | 110+10 mm |
| XL | >110 mm |



- | |
|-------------|
| min. 240 mm |
| min. 240 mm |
| min. 240 mm |
| min. 240 mm |



UK & Europe HQ
Supermax Healthcare Limited
Units 12-16 Titan Drive,
Fengate, Peterborough,
PE1 5XN, UK

+441733 615215
sales@supermax.co.uk
www.aureliagloves.co.uk



Europe
Supermax Healthcare (Europe) Limited
Units 1001-1012
Building 1000, Gateway Business Park,
New Mallow Road Cork, T23 F2H7, Ireland

+353 15256724
sales@supermax.ie
www.aureliagloves.ie




@aureliaglovesuk