



AURELIA® UNIQUE TPE BLUE

Powder Free TPE Examination Gloves

AURELIA® GLOVES
A SUPERMAX HEALTHCARE BRAND

sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building 1000,
Gateway Business Park,
New Mallow Road, Cork
Ireland

aureliagloves.ie



Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 455 1-4 Medical gloves for single use. BS EN 1186:2002 Materials and articles in contact with foodstuffs.



10 x 200

S M L XL

48226 48227 48228 48229

Sizes & Dimensions

Single wall thickness

Physical Properties

S 103+2 mm	min. 270 mm
M 105+2 mm	min. 270 mm
L 115+2 mm	min. 270 mm
XL > 120 mm	min. 270 mm



- ① min. 0.06 mm
- ② min. 0.06 mm
- ③ min. 0.05 mm

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 400%	Min. 400%
Strength at Break	Min. 3.0N	Min. 3.0N



TPE



Powder Free



3 Years



Smooth



2.2g +/- 0.2g



Blue

