



# AURELIA® DELIGHT CLEAR PF

Premium Powder Free Vinyl Examination Gloves

## AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building 1000,  
Gateway Business Park,  
New Mallow Road, Cork,  
Ireland

[aureliagloves.ie](http://aureliagloves.ie)



Regulations	PPE Regulation (EU) 2016/425
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 455 1-4 Medical gloves for single use. BS EN 1186:2002 Materials and articles in contact with foodstuffs.





S	M	L	XL
38226	38227	38228	38229

### Sizes & Dimensions

### Single wall thickness

### Physical Properties

	
S 80+10 mm	min. 240 mm
M 95+10 mm	min. 240 mm
L 110+10 mm	min. 240 mm
XL > 110 mm	min. 240 mm



- ① min. 0.05 mm
- ② min. 0.05 mm
- ③ min. 0.04 mm

	Unaged	Aged
Tensile Strength	Min. 11.0 MPa	Min. 11.0 MPa
Elongation	Min. 350%	Min. 350%
Strength at Break	Min. 3.6N	Min. 3.6N



Vinyl



Powder Free



3 Years



Smooth



4.5g +/- 0.2g



Clear

