



# AURELIA® UNIQUE TPE CLEAR

Powder Free TPE Examination Gloves

## AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building 1000,  
Gateway Business Park,  
New Mallow Road, Cork,  
Ireland

[aureliagloves.ie](http://aureliagloves.ie)



Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 455 1-4 Medical gloves for single use. BS EN 1186:2002 Materials and articles in contact with foodstuffs.



Sizes & Dimensions		Single wall thickness	Physical Properties													
 S 103+2 mm M 105+2 mm L 115+2 mm XL > 120 mm	 min. 270 mm min. 270 mm min. 270 mm min. 270 mm	 ① min. 0.06 mm ② min. 0.06 mm ③ min. 0.05 mm	<table border="1"> <thead> <tr> <th></th> <th>Unaged</th> <th>Aged</th> </tr> </thead> <tbody> <tr> <td>Tensile Strength</td> <td>Min. 14.0 MPa</td> <td>Min. 14.0 MPa</td> </tr> <tr> <td>Elongation</td> <td>Min. 400%</td> <td>Min. 400%</td> </tr> <tr> <td>Strength at Break</td> <td>Min. 3.0N</td> <td>Min. 3.0N</td> </tr> </tbody> </table>		Unaged	Aged	Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa	Elongation	Min. 400%	Min. 400%	Strength at Break	Min. 3.0N	Min. 3.0N	
	Unaged	Aged														
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa														
Elongation	Min. 400%	Min. 400%														
Strength at Break	Min. 3.0N	Min. 3.0N														



TPE



Powder Free



3 Years



Smooth



2.2g +/- 0.2g



Clear

