



AURELIA® IGNITE

7mil, Nitrile Powder Free Examination Gloves

AURELIA® GLOVES

A SUPERMAX HEALTHCARE BRAND

sales@supermax.ie

+353 15256724

Units 1001 – 1012, Building 1000,
Gateway Business Park,
New Mallow Road, Cork,
Ireland

aureliagloves.ie

Regulations	PPE Regulation (EU) 2016/425 Medical Device Regulation (EU) 2017/745, 93/42/EEC
Quality Standards	ISO 9001:2015 Quality Management Systems ISO 13485: 2016 Medical Devices Quality Management Systems
Harmonised Standards	BS EN ISO 374 - 1, 2, 4, 5 Protective gloves against dangerous chemicals and microorganisms. BS EN 455 1-4 Medical gloves for single use. EN 420: 2003 + A1: 2009 Protective gloves. General requirements and test methods. BS EN 1186: 2002. Materials and articles in contact with foodstuffs.

10 x 100
 * 10 x 90

M L *XL *XXL
 97887 97888 97889 97880

Sizes & Dimensions	Single wall thickness	Physical Properties
--------------------	-----------------------	---------------------

M	95+10 mm	min. 240 mm
L	110+10 mm	min. 240 mm
XL	120+10 mm	min. 240 mm
XXL	> 120 mm	min. 240 mm

①	min. 0.25 mm
②	min. 0.20 mm
③	min. 0.09 mm

	Unaged	Aged
Tensile Strength	Min. 14.0 MPa	Min. 14.0 MPa
Elongation	Min. 500%	Min. 400%
Strength at Break	Min. 6N	Min. 6N

Nitrile	Powder Free	5 Years
Diamond	8.4g +/- 0.2g	Orange

